

SAFEGUARD MEDICAL™

PIN®



Overview

The PIN® provides rapid intraosseous (IO) access when intravenous (IV) access is challenging or impossible. The PIN® can be used in emergency medical or trauma cases and can remain in-situ for up to 24 hours. The PIN® can be used for IO access in the humeral head, proximal and distal tibia sites when IV access is difficult or impossible to obtain in emergent, urgent, or medically necessary cases for up to 24 hours.

Key Benefits

- Rapid and easy IO insertion.
- Single needle solution - one size for all, including paediatrics.
- 15-Gauge needle, pharmacokinetics of which are well established.
- Stabilisation accessory supplied.
- Single-use product - no driver cleaning required between uses.
- No power or battery required for insertion.
- Lightweight and robust
- Indicated for access to six IO sites
- Latex-free and non-pyrogenic
- Product code: PDF103

Specifications



Weight: 4.8kg

For further information on the PIN®, visit: www.safeguardmedical.com